Sibelmed

LOW COST

CABINA S-BASIC AUDIOMETRIC TEST BOOTH



- Great noise attenuation, 41 db SPL at 10KHz
- Quick and easy assembly
- Compatible with the main audiometers of the market
- Connection board included







S-BASIC · audiometric test booth

Description

The audiometric test booths **S-BASIC** has been designed by SIBEL S.A.'s RDI department using the newest materials in noise attenuation. It has allowed to reduce its weight, get a high isolation index, great functionality and pleasant looking.

The external surface is made of washable pannels and the interior is pleasantly illuminated and comfortable.

Features

External booth: 96 x 96 x 208 cm.Internal booth: 86 x 86 x 207 cm.

• Door measures: 70 x 180 cm.

Window: 54 x 46 cm.
Foldable table: 80 x 40 cm.
External material: Melamine
Internal material: Pyramid foam
Interior illumination: 35W lamp
Power supply: 120-220V at 50-60 Hz

Connection boards: 5+1Weight: <200 kg

* Optional

TWO DIGITAL CHANNELS AUDIOMETER

Suitable to your needs: 6 models*.

- Advanced digital technology based on DSP.
- Calculation of hearing loss and diagnosis.
- Internal database > 1.000 tests.
- Suprathreshold tests.
- High frequency.
- Musical frequencies.
- Automatic masking.
- Connectivity with other management systems.
- · Customizable.
- Audiometry sofware W50

See configurations*

Regulations


The Audiometric Test Booth S-BASIC complies with:

• EN 60601-1:2006 + AC:2010, EN 60601-1-6:2010, EN 62366:2008, EN ISO 14971:2012, EN 1041:2008, EN ISO 15223-1:2012

This product has been validated in laboratory according to UNE-EN ISO 11957:1997 regulation and calculation of the global index of acoustic reduction in according with UNE EN ISO 717-1:1997.

The process of designing, production and customer service will be made according to ISO 9001:2008, EN 13485:2003 + AC:2007 quality standards and to the 93/42/CEE european directive of Sanitary Products requirements.

Noise attenuation





SIBEL S.A., Rosselló 500, 08026 Barcelona (SPAIN)

National Sales: Tel. +34 93 436 00 08 e-mail: comercial@sibelgroup.com Export Sales: Tel.+34 93 436 00 07 e-mail: export@sibelgroup.com

Technical support: +34 93 433 54 50 e-mail: sat@sibelgroup.com Fax: +34 93 436 16 11

SIBEL S.A. belongs to Sibel GROUP

www.sibelmed.com • www.sibelgroup.com







