



MANUFACTURER N°1

PULMONARY FUNCTION MOUTHPIECES



- Class IIa certified according to 93/42/CEE Directive.
- Different measures for any spirometer in the market.

28 ext. x 60mm 30 ext. x 60mm 30,5 ext. x 65mm 22 ext. x 60mm
 27 ext. x 70mm 21 ext. x 70mm 26,0 ext. x 73mm ...

- Available in two presentation formats (100u bags.)



In bulk (Standard)



Single packaging (Optional)

CODE	MEASURES	USE	COMPATIBILITY
01555	28 ext. x 60mm	ADULT	SIBELMED (Datospir), LETI (Peak Flow PF-Control)
01554	30 ext. x 60mm	ADULT	SIBELMED (Datospir, PEAK 10), CLEMENT CLARKE, MICROMEDICAL (Microlab), MIR (Spirolab, Spirobank, ...)
01553	30,5 ext. x 60mm	ADULT	COSMED (Pony), VITALOGRAPH
01552	27 ext. x 70mm	ADULT	POCH-MILLAS, MINATO (Autospiro AS500), EYMASA (Neumos-12)
01557	21 ext. x 70mm	PEDIATRIC	SIBELMED (Datospir), POCH-MILLAS
01551	26 ext. x 73mm	ADULT	FUKUDA, JAEGER, WELCH-ALLYN
01558	22 ext. x 60mm	PEDIATRIC	COSMED (Pony,...), CLEMENT CLARKE, VITALOGRAPH, MICRO MEDICAL (Co-oximeters MICRO-CO, Smoke Check)

FOR YOUR SAFETY ALWAYS REQUIRE A CERTIFICATE GIVEN BY AN AUTHORIZED ORGANIZATION; DO NOT ACCEPT A SIMPLE AUTOCERTIFICATION.



NB Identification N° of the Spanish Drugs and Sanitary Products Agency.

OTHER ACCESSORIES, SPARE PARTS & CONSUMABLES



Peak flow



Disposable transducers



Calibration syringe



Bacterial filter



Nose clip



Paper