

# Hemo Control Manager

## Hemoglobin & Hematocrit Analyzer



Lab accurate hemoglobin results with data management



(Printer not included)

Hemo Control Manager has the same standard features as Hemo Control

Barcode identification of patient ID, operator ID and QC

Connectivity to PCs via cable or Bluetooth® technology

Memory for 4,000 patients results and 500 QC results

Quality control reminders

Flags results out of range or rejected results

Access control and QC lock-out function

### Connected hemoglobin testing

The Hemo Control analyzer features a bi-directional interface using the LIS2-A2 communication standard.

### Ordering information

#### Pack includes

Analyzer, barcode scanner, power adaptor, user guide, quick guide, optic cleaner, control cuvette and wall chart.

#### Connectivity options

**Hemo Control** can be upgraded with Data Management functions using the Add Pack Hemo Control DM option.

**Hemo Control Manager** supports Data Management functions at delivery and features integrated Bluetooth connectivity.

**Optional PC software** If you are not planning to connect to a POC Data Manager/LIS directly or just require a PC connection, the Hemo Connect software allows users to configure data management functions and download test results when the test is performed or later from the device memory.

Description	Reference No.	Contents
Hemo Control Manager	3041-0010-0229	1 pack
Hemoglobin microcuvette (canisters)	3000-3012-0765	1 canister of 50 pcs (4 x 50 pcs per pack)
Hemoglobin microcuvette (Individually packed)	3000-3013-0278	1 box of 50 pcs
Control solution Hb-con low	3000-6121	1 pc
Control solution Hb-con set 2	3000-6128	1 set
Bluetooth dongle	3041-7211-0505	1 pc
Bluetooth printer MCP 1880	3040-7011-0452	1 pc
Cleaner	3000-6232	1 pack of 5 cleaners
HemoConnect Software	3040-7231-0132	1 pack