

Technical data:

MPD 250K100M Professional Step-On personal floor scale with approval for medical use

Volume / materials / weights

Dimensions of display device (WxDxH):	200 mm x 55 mm x 128 mm
Dimensions of weighing plate (WxDxH):	365 mm x 360 mm x 80 mm
Gross weight:	10,500 kg
Net weight:	10 kg
Overall dimensions mounted (WxDxH):	365 mm x 490 mm x 120 mm

Service

DKD Certificate:	963-129
Verification:	965-129

Technical data about verification

n/a:	200000000
Verification approval conforming to 90/384/EEC:	ja
Verification by KERN possible:	ja
Verification class:	III
Verification value:	100 g

Balance features

Display digit height:	2,500 cm
Level indicator:	ja
Revolving screw feet:	ja
Weighing units:	kg

Category

Category:	Balances
-----------	----------

Recommended adjusting weight:	200 kg 0
-------------------------------	----------

Measuring system

Linearity:	100 g
Minimum load:	2 kg
Readout [d]:	100 g
Reproducibility:	100 g
Weighing range [Max]:	250 kg
Weighing system:	Strain gauge

Power supply	
Charging time:	12 h
Input voltage:	100 V - 240 V 50
Operating time (Backlight off):	40 h
Operating time (Backlight on):	20 h
Environmental conditions	
n/a:	80 %
n/a:	40 °C
n/a:	0 °C
Packaging & shipping	
Dimensions packaging (WxDxH):	400 mm x 510 mm x 100 mm